

# Ultra

The 8-place Ultra centrifuge has a unique, space-saving triangle shape to fit easily into a corner. This unit, made of sturdy aluminum alloy, accommodates most test tubes from 3ml-15ml. The wide stance eliminates wobble, while the suction-cupped feet prevent slipping. The brushless motor (110V only) is virtually maintenance-free and the clear lid allows easy viewing of specimens while in use.

Choose from 2 different models to meet your laboratory separation needs.

## ADVANTAGES

- Accurate speeds for specific g-forces
- Proper separations without cell damage
- Ideal for blood, urine and other fluids
- Quiet, vibration-free operation
- Unique space-saving triangle design

## FEATURES

- 8-tube capacity - 3ml to 15ml size tubes
- Auto-off 30-minute timer with bell
- Brushless motor (110V only) - maintenance free
- Blue LED power indicator light
- Suction-cupped feet prevent slipping
- All-metal one-piece body
- 110V power (220V available)

## ULTRA SELECT

Selectable 4-position Speed Control

- 3,300 rpm for blood - 1,350 g
- 1800 rpm for urines and semen - 400 g
- 1,500 rpm for fecals - 275 g
- 800 rpm for other tests - 80 g

Timer and lid safety shut-off switch

## ULTRA VARIABLE

3,300 rpm variable speed - blood, urines, etc.  
Timer and lid safety shut-off switch  
Maximum RCF=1,350 g

## INCLUDES

- Eight-place tube rotor
- Eight 15ml tube sleeves (TS-15ml)
- Eight 13x75mm tube sleeve inserts (TS-1375)



## Ultra

Model #	Description
ULC-08AV-1501	Variable Speed
ULC-08AS-1501	4-speed Selectable

- All-Metal One-Piece Design
- 2 Models to Meet Any Lab Need
- Locking Lid for Safety

## WEIGHT AND DIMS

Height: 9" (230mm)  
Length: 13" (330 mm)  
Width: 13" (330mm)  
Weight: 15.5 lbs (7kg)



P 770.270.1394  
F 770.270.2389

865 Marathon Parkway  
Lawrenceville GA 30046

LWScientific.com