

BaXstrap® Spineboard

As a rescuer you want your spineboard to be comfortable for our patient, rigid and lightweight. Further you also want to store the board easily. The BaXstrap Spineboard provides this and is one of the lightest boards on the market. The contoured design allows the board to nest compactly. This unique feature also offers patients more comfort, making it less likely that they would move around and aggravate their injuries.



- Dimensions: 16" x 72" x 2.5"
- Weighs 13.5 lbs. (single board)
- Thickness of 2.3 inches
- 14 handholds
- Handhold width of 2 inches
- Handhold length of 5.25 inches
- Minimum storage space for 2 boards of approximately 4.00 inches
- 12 Pins
- Carbon fiber tubes
- The spineboard shall be CE marked
- Operating temperature: -18°C to 43°C (0°F to 110°F)
- Storage temperature: -18°C to 43°C (0°F to 110°F)
- Limited lifetime warranty

- No seams make it easy to clean and decontaminate
- Extra large handholds make it easy to grab, even with heavy gloved hands
- Raised handholds make it easy to pick up
- Pins are molded into the board, maximizing their strength
- Includes strap holes for pediatric patients, eliminating the need to use towels to fill the extra space between patient and strap
- Can be used for water rescue, in the mountains, at home or on the road
- Durable molded-in customized labeling available
- Radiolucent for clear view of the spine, MRI and CT scan compatible
- Latex-free