

## Laerdal Silicone Masks



**Your patients come in different shapes and sizes and each has their own unique facial anatomy. Effective mask seal depends on the ability of your mask to be flexible. Laerdal Silicone Masks give you that necessary flexibility.**

- Silicone masks have a one-piece design.
- Soft silicone cuff and dome allows for airtight seal while providing patient comfort.
- Transparent for high visibility of mouth and nose.
- Fits standard 15/22mm connectors.
- Multi-Function Mask Cover allows for versatile use of Adult and Child Masks.
- Reusable and latex-free