

Laerdal Silicone Resuscitators (LSR)

Over the past 35 years, Laerdal Silicone Resuscitators have proven to be a leader in reusable silicone resuscitators. The new design enhances convenience and maintains the outstanding quality, efficiency and economy that healthcare professionals have come to rely on.



Description of the LSR features:

- Transparent silicone material provides excellent bag re-expansion and is resistant to extreme temperatures
- Superior feel of lung compliance helps reduce risk of over inflation
- Available in Adult, Pediatric and Preterm sizes.
- Integrated Intake/Reservoir Valve for efficient oxygen delivery and ease of cleaning.
- See through one-piece Laerdal Silicone Mask with Multi-Function Mask Cover Material allows for various cleaning, disinfecting and sterilization methods.
- Pressure Relief Valve setting on Pediatric and Preterm models.
- Latex-free components.
- Various configuration and storage container options.
- Private label and storage container options available.