

STAINLESS STEEL CLOSED CASE CARTS

Case carts are the ultimate choice for sterility control during transportation.

COMMON FEATURES:


STAINLESS STEEL CLOSED CASE CARTS

- ✓ Sealed doors safeguard contents against cross-contamination
- ✓ Positive door latches secure doors shut to keep contents safe and sterile
- ✓ Hold open latch secures door when open for easy cleaning
- ✓ Stainless steel is easy to clean and sanitize, resists corrosion and stains
- ✓ All-welded construction eliminates nuts and bolts for a stronger structure and easier-to-clean surface
- ✓ Extended push handle offers more control and increases interior space
- ✓ Locking door latch available for added security




#6950

Shown with optional locking door latch



Floor is cross-broke and angled to middle with small holes to drain



Pull-out shelves equipped with locking mechanism



Door latch



CLOSED CASE CARTS

MODEL	ACTUAL SHELF WIDTH	EXTERNAL DIMENSIONS			INTERNAL DIMENSIONS			SHIP WEIGHT (LBS.)
		W	D	H	W	D	H	
ONE DOOR								
6930	24"	36"	29"	39"	25 1/2"	25"	29 1/4"	229
6932	24"	36"	29"	54"	25 1/2"	25"	44 1/4"	250
6940	18"	30"	29"	39"	19 1/2"	25"	29 1/4"	197
6942	18"	30"	29"	54"	19 1/2"	25"	44 1/4"	220
TWO DOOR								
6947†	31 3/8"	43"	29"	54"	34"	25"	44 1/4"	291
6952†	31-3/8"	43"	29"	39"	34"	25"	29 1/4"	257
6935†	36"	48"	29"	39"	38 1/4"	25"	29 1/4"	282
6937†	36"	48"	29"	54"	38 1/4"	25"	44 1/4"	309
6945	13 1/2"	43"	29"	54"	15 1/4"*	25"	44 1/4"	296
6950	13 1/2"	43"	29"	39"	15 1/4"*	25"	29 1/4"	262

† Denotes no center post *Two sections/center post



STAINLESS STEEL CLOSED CASE CARTS (CONTINUED)



#6935



Model 6935 shown with one door secured in open position with hold open latch



Model 6935 with both doors open

CLOSED CASE CART SPECIFICATIONS

CART	Non-magnetic stainless steel
FINISH	#3 Finish
MATERIAL	16-gauge stainless steel
BUMPER	Replaceable, stand-out corner bumper, non-marring
SHELVES	Each model includes one stainless steel wire pull-out shelf, except models 6945 and 6950, which include 2 pull-out shelves; adjustable on 2" centers
CASTERS	5" stainless steel casters are maintenance-free



OPTIONAL, ADDITIONAL FIXED WIRE SHELVES

MODEL	DESCRIPTION	COMPATIBILITY
3810	13 1/2", 75 lb capacity	Fits 6945 and 6950
3820	18", 75 lb capacity	Fits 6940 and 6942
3830	24", 125 lb capacity	Fits 6930 and 6932
3840	31 3/8", 125 lb capacity	Fits 6947 and 6952
3850	36", 125 lb capacity	Fits 6935 and 6937



OPTIONAL, ADDITIONAL PULL-OUT WIRE SHELVES

MODEL	DESCRIPTION	COMPATIBILITY
3814	13 1/4", 75 lb capacity	Fits 6945 and 6950
3824	18", 75 lb capacity	Fits 6940 and 6942
3834	24", 125 lb capacity	Fits 6930 and 6932
3844	31 1/2", 125 lb capacity	Fits 6947 and 6952
3854	36", 125 lb capacity	Fits 6935 and 6937



STAINLESS STEEL OPEN CASE CARTS

Open case carts can be configured to meet your every need!

- Maximum design flexibility for your best fit
- Multiple uses as case carts, supply carts, distribution carts, linen carts, etc.
- Stainless steel frame is easy to clean and sanitize
- Solid stainless steel or wire shelves allow options for different size supplies

- All-welded construction eliminates nuts and bolts for a stronger structure and easier-to-clean surface
- Fits elevators and adaptable for cart pick-up systems
- 6" Casters for easy maneuverability
- Easy to configure with 5 1/2" (140) space between shelf rungs
- No tools required

BASE MODEL 7910

(38 1/2"W x 28 3/4"D x 38 1/2"H) Accommodates up to 5 shelves. Select all wire, all stainless or a combination of wire and stainless.

7910W	Wire Shelf
7910SS	Stainless Steel Shelf
7910SD	Solid Shelf with Divider (3 adjustable dividers with hardware)

BASE MODEL 7920

(42 1/2"W x 30 3/4"D x 44"H) Accommodates up to 6 shelves. Select all wire, all stainless or a combination of wire and stainless.

7920W	Wire Shelf
7920SS	Stainless Steel Shelf
7920SD	Solid Shelf with Divider (3 adjustable dividers with hardware)

BASE MODEL 7930

(32 1/2"W x 32 3/4"D x 61 1/2"H) Accommodates up to 9 shelves. Select all wire, all stainless or a combination of wire and stainless.

7930W	Wire Shelf
7930SS	Stainless Steel Shelf
7930SD	Solid Shelf with Divider (3 adjustable dividers with hardware)

BASE MODEL 7940

(52 1/2"W x 27 1/2"D x 61 1/2"H) Accommodates up to 9 shelves. Select all wire, all stainless or a combination of wire and stainless.

7940W	Wire Shelf
7940SS	Stainless Steel Shelf
7940SD	Solid Shelf with Divider (3 adjustable dividers with hardware)



#7910



#7920
shown with optional bins
(120118 gray, 120119 red, 120120 blue)



#7930
shown with optional bins
(120120 blue, 120119 red)



#7940
shown with optional bins
(120119 red, 120118 gray, 120120 blue)



OPEN CASE CART ACCESSORIES

Choose from shelf dividers, tote box drawers, tote box drawer dividers, stainless steel casters, Staph-Chek® vinyl covers, and waste bag attachment.

STANDARD PRODUCTS
SHIP IN **10** DAYS OR LESS

CART COVERS

Navy is standard. Also available in Light Blue, Mint Green, Mauve, Tan and White.

29650	Use with Base Model 7910
29652	Use with Base Model 7920
29655	Use with Base Model 7930
29660	Use with Base Model 7940



#7930
shown with optional 29655 cart cover

PLASTIC BINS

Available in Red, Blue or Gray. Please indicate your bin color choice at time of order.

120118	Small Bins (8 1/4"W x 10 7/8"D x 5"H)
120119	Medium Bins (10 7/8"W x 16 1/2"D x 6"H)
120120	Large Bins (17 3/8"W x 22 3/4"D x 6"H)



#120120 blue bin

BIN SLIDE BRACKETS

120110	Bin Bracket Kit (pair) Use with Base Models 7910, 7920, 7940
120111	Bin Bracket Kit (pair) Use with Base Model 7930

BAG KITS

120112	Trash Bag Kit (Bag and Bracket) Use with Model 7910
120113	Trash Bag Kit (Bag and Bracket) Use with Models 7920, 7940
120114	Trash Bag Kit (Bag and Bracket) Use with Model 7930
53511	Replacement Bag, Blue Staph-Chek

STAINLESS STEEL CASTERS

120125	5" Stem Caster
120126	5" Stem Caster with Brake
120128	6" Stem Caster
120127	6" Stem Caster with Brake



#7920
shown with optional 120113 trash bag kit