




*Varied sizes to stabilize the fingers,
wrists and legs*

ORT32400

Aluminum Finger Splints

- A** Medline's Padded Finger Strips can be cut to size
- Malleable aluminum lined with white hydroloc foam

Description	Item Number	Qty
½" (1.3 cm) L x 9" (22.9 cm) W	ORT3210019	12/cs
¾" (1.9 cm) L x 9" (22.9 cm) W	ORT3210029	12/cs
½" (1.3 cm) L x 18" (45.7 cm) W	ORT32100118	12/cs
¾" (1.9 cm) L x 18" (45.7 cm) W	ORT32100218	12/cs
1" (2.54 cm) L x 9" (22.9 cm) W	ORT3210039	12/cs
1" (2.54 cm) L x 18" (45.7 cm) W	ORT32100318	12/cs
1½" (3.8 cm) L x 18" (45.7 cm) W	ORT32100418	12/cs



A



B

Gutter/ Spoon/ Curved Finger Splint protects the end of the finger

- Malleable aluminum lined with blue polyester urethane foam
- Curved to protect the end of the finger from further injury

Description	Item Number	Qty
1½" (3.81 cm) L x ⅝" (2.3 cm) W, Small	ORT32300S	12/cs
3" (7.62 cm) L x ⅝" (2.3 cm) W, Medium	ORT32300M	12/cs
6" (15.2 cm) L x ⅝" (2.3 cm) W, Large	ORT32300L	12/cs



C

- C** Medline's Frog Finger Splint folds around finger to stabilize
- Malleable aluminum lined with blue polyester urethane foam

Description	Item Number	Qty
2¼" (5.72 cm) L x 2¾" (7 cm) W, Small	ORT32200S	12/cs
2½" (6.35 cm) L x 3" (7.62 cm) W, Medium	ORT32200M	12/cs
3" (7.62 cm) L x 3½" (8.9 cm) W, Large	ORT32200L	12/cs



D

- D** Medline's Four-Prong Finger Splint fully encloses the injured finger

- Malleable aluminum lined with blue polyester urethane foam

Description	Item Number	Qty
1½" (3.81 cm) L, Small	ORT32600S	12/cs
2¾" (7 cm) L, Medium	ORT32600M	12/cs
3¼" (8.26 cm) L, Large	ORT32600L	12/cs

* Measurements based upon dimensions similar to those of patient's injured finger

Aluminum Finger Splints



E Fold-Over Finger Cot protects the fingertip and finger

- Malleable aluminum lined with blue polyester urethane foam

Description	Item Number	Qty
1½" (3.81 cm) L, X-Small	ORT32400XS	12/cs
2¼" (5.72 cm) L, Small	ORT32400S	12/cs
3¼" (8.26 cm) L, Medium	ORT32400M	12/cs
3¾" (9.35cm) L, Large	ORT32400L	12/cs



G Medline's Baseball Finger Splint stabilizes fractures of the distal phalanx

- Malleable aluminum lined with blue polyester urethane foam

Description	Item Number	Qty
3" (7.6 cm) L, Small	ORT32500S	12/cs
4¼" (10.8 cm) L, Medium	ORT32500M	12/cs
5" (12.7 cm) L, Large	ORT32500L	12/cs



F Assorted Finger Splint Kit

- An assortment of 5 styles of padded finger splints including frog, spoon, cot, 4-prong and baseball
- Contains 3 of each style and size
- Packaged in a convenient, clear plastic storage case

Description	Item Number	Qty
Finger Splint Kit	ORT32000	ea.

* Measurements based upon dimensions similar to those of patient's injured finger

Finger Splints




Questions? Comments?
Call 1-800-MEDLINE, visit our Web site, www.medline.com, or contact your Medline sales representative.



A Medline's Stackies Kit and Replacements offer quick therapy

- Molded polypropylene splints protect and stabilize the distal fingertip
- Kit contains an assortment of 30 splints packaged in clear plastic storage container
- Available in eight sizes
- Open nail design with ventilation holes

Description	Item Number	Qty
Stackies Kit, Assorted Sizes	ORT32700	ea.
Replacement, Size 1	ORT327001C	6/cs
Replacement, Size 2	ORT327002C	6/cs
Replacement, Size 3	ORT327003C	6/cs
Replacement, Size 4	ORT327004C	6/cs
Replacement, Size 5	ORT327005C	6/cs
Replacement, Size 5.5	ORT3270055C	6/cs
Replacement, Size 6	ORT327006C	6/cs
Replacement, Size 7	ORT327007C	6/cs



B Plastalume Finger Splints

- Kit contains six of each size splint
- Aluminum splint is lined with foam and covered in white vinyl
- Bulbs on some styles reduce the possibility of the splint turning or rolling
- Radiolucent
- Measurements based upon dimensions similar to those of patient's injured finger

Description	Item Number	Qty
Finger Splint Kit, Assorted Sizes	ORT328000	ea.
Finger Splint, 1¼" (3.2 cm) L	ORT328001	12/cs
Finger Splint, 2¼" (5.72 cm) L	ORT328002	12/cs
Finger Splint, 3¼" (8.26 cm) L	ORT328003	12/cs
Finger Splint, 4¼" (10.8 cm) L	ORT328004	12/cs
Finger Splint w/ Bulb, 3¼" (8.26 cm) L	ORT328005	12/cs
Finger Splint w/ Bulb, 4¼" (10.8 cm) L	ORT328006	12/cs
Finger Splint w/ Bulb, 5¼" (13.3 cm) L	ORT328007	12/cs
Finger Splint w/ Bulb, 6¼" (15.9 cm) L	ORT328008	12/cs
Finger Splint w/ Bulb, 7¼" (18.4 cm) L	ORT328009	12/cs



C Medline's Aluminum Colles Splint

- Padded with blue polyester urethane foam for comfort
- Pre-molded for an anatomical fit
- Malleable to allow exact positioning



D Medline's Aluminum Forearm Splint

- Padded with blue polyester urethane foam for comfort
- Pre-molded for an anatomical fit
- Malleable to allow exact positioning

Description	Item Number	Qty
Right		
7½" (19.1 cm) L, 2¾" (7 cm) W, Small	ORT33100RS	ea.
9½" (24.1 cm) L, 3¾" (9.5 cm) W, Medium	ORT33100RM	ea.
10½" (26.7 cm) L, 4¼" (10.8 cm) W, Large	ORT33100RL	ea.
Left		
7½" (19.1 cm) L, 2¾" (7 cm) W, Small	ORT33100LS	ea.
9½" (24.1 cm) L, 3¾" (9.5 cm) W, Medium	ORT33100LM	ea.
10½" (26.7 cm) L, 4¼" (10.8 cm) W, Large	ORT33100LL	ea.

Description	Item Number	Qty
Right		
7" (17.8 cm) L, 3½" (8.9 cm) W, Small	ORT33300RS	ea.
8½" (21.6 cm) L, 3½" (8.9 cm) W, Medium	ORT33300RM	ea.
9" (22.9 cm) L, 4½" (11.4 cm) W, Large	ORT33300RL	ea.
Left		
7" (17.8 cm) L, 3½" (8.9 cm) W, Small	ORT33300LS	ea.
8½" (21.6 cm) L, 3½" (8.9 cm) W, Medium	ORT33300LM	ea.
9" (22.9 cm) L, 4¼" (10.8 cm) W, Large	ORT33300LL	ea.

* Measurements based upon dimensions similar to those of a particular patient's injured limb

Aluminum Splints



A Medline's Aluminum Metacarpal Splint

- Padded with blue polyester urethane foam for comfort
- Pre-molded for an anatomical fit
- Malleable to allow exact positioning

Description	Item Number	Qty
Right		
8" (20.3 cm) L, 3" (7.62 cm) W, Small	ORT33200RS	ea.
10½" (26.7 cm) L, 3¾" (9.5 cm) W, Medium	ORT33200RM	ea.
11½" (29.2 cm) L, 4½" (11.4 cm) W, Large	ORT33200RL	ea.
Left		
8" (20.3 cm) L, 3" (7.62 cm) W, Small	ORT33200LS	ea.
10½" (26.7 cm) L, 3¾" (9.5 cm) W, Medium	ORT33200LM	ea.
11½" (29.2 cm) L, 4½" (11.4 cm) W, Large	ORT33200LL	ea.



B Medline's Aluminum Cock-Up Splint

- Padded with blue polyester urethane foam for comfort
- Pre-molded for an anatomical fit
- Malleable to allow exact positioning

Description	Item Number	Qty
6½" (16.5 cm) L, 2½" (6.35 cm) W, Small	ORT33500S	ea.
8" (20.3 cm) L, 3" (7.62 cm) W, Medium	ORT33500M	ea.
9" (22.9 cm) L, 3" (7.62 cm) W, Large	ORT33500L	ea.

* Measurements based upon dimensions similar to those of a particular patient's injured limb



C

C Medline's Aluminum Mason Allen Splint

- Padded with blue polyester urethane foam for comfort
- Pre-molded for an anatomical fit
- Malleable to allow exact positioning

Description	Item Number	Qty
9½" (24.1 cm) L, 5" (12.7 cm) W, Small	ORT33600S	ea.
12½" (31.8 cm) L, 6" (15.2cm) W, Large	ORT33600L	ea.



D

D Medline's Posterior Tibia and Fibula Aluminum Splints

- Available in knee or femoral length
- Padded for patient comfort
- Fits right or left leg

Description	Item Number	Qty
Knee Length		
11" (27.9 cm) L, 7" (17.8 cm) W, Small	ORT34900S	ea.
13½" (34.5 cm) L, 8" (20.3 cm) W, Medium	ORT34900M	ea.
16" (40.6 cm) L, 9½" (24.1 cm) W, Large	ORT34900L	ea.
Femoral Length		
16" (40.6 cm) L, 5" (12.7 cm) W, Small	ORT34800S	ea.
20" (50.8 cm) L, 6" (15.2 cm) W, Medium	ORT34800M	ea.
23" (58.4 cm) L, 6½" (16.5 cm) W, Large	ORT34800L	ea.

* Measurements based upon dimensions similar to those of a particular patient's injured limb