

Cast protection and foot support for any weather

ORT30300

Cast Boots



A) Black Canvas Cast Sandal offers a lightweight way to protect your cast

- Made of a quality canvas duck upper for durability
- Two hook-and-loop closures over the top of the foot and one above the heel for a secure fit
- Thick rocker sole is made of a lightweight EVA
- Open toe and open heel for comfort and ventilation
- Heel and toe pitch promotes a normal gait
- XXS size comes in pediatric print

Description	Item Number	Qty
5½" (14 cm) L, 3" (7.6 cm) W, 2X-Small Shoe size: Child's 5½-7½	ORT29100XXS	ea.
7¼" (18.4 cm) L, 3½" (8.9 cm) W, X-Small Shoe size: Women's 3-4½, Child's 10½-12½	ORT29100XS	ea.
8¾" (22.2 cm) L, 4½" (11.4 cm) W, Small Shoe size: Women's 5-7	ORT29100S	ea.
10½" (26.7 cm) L, 4¾" (12.1 cm) W, Medium Shoe size: Men's 9-10, Women's 10-12	ORT29100M	ea.
12" (30.5 cm) L, 5¼" (13.3 cm) W, Large Shoe size: Men's 12-15½	ORT29100L	ea.
13½" (34.3 cm) L, 5½" (14 cm) W, X-Large Shoe size: Men's 16+	ORT29100XL	ea.



B) Vinyl Closed Toe Cast Boot protects cast from cold, damp weather

- Closed toe and heel is ideal for colder climates
- Vinyl upper resists water and easily wipes clean after exposure to the elements
- Two hook-and-loop closures over the top of the foot and one above the heel for a secure fit
- Thick rocker sole is made of a lightweight EVA
- Heel and toe pitch promotes a normal gait

Description	Item Number	Qty
5½" (14 cm) L, 3" (7.6 cm) W, 2X-Small Shoe size: Child's 5½-7½	ORT29210XXS	ea.
7¼" (18.4 cm) L, 3⅝" (9.2 cm) W, X-Small Shoe size: Women's 3-4½, Child's 10½-12½	ORT29210XS	ea.
8¾" (22.2 cm) L, 4¾" (11.1 cm) W, Small Shoe size: Women's 5-7	ORT29210S	ea.
10½" (26.7 cm) L, 4¾" (12.1 cm) W, Medium Shoe size: Men's 9-10, Women's 10-12	ORT29210M	ea.
12" (30.5 cm) L, 5¼" (13.3 cm) W, Large Shoe size: Men's 12-15½	ORT29210L	ea.
13½" (34.3 cm) L, 5½" (14 cm) W, X-Large Shoe size: Men's 16+	ORT29210XL	ea.

C) Vinyl Open Toe Cast Boot protects cast while providing some ventilation

- Closed heel design helps keep the cast dry in damp weather
- Open toe for comfort and breathability
- Vinyl upper resists water and easily wipes clean after exposure to the elements
- Two hook-and-loop closures over the top of the foot and one above the heel for a secure fit
- Thick rocker sole is made of a lightweight EVA
- Heel and toe pitch promotes a normal gait

Description	Item Number	Qty
4½" (11.4 cm) L, 3½" (8.9 cm) W, 2X-Small, Shoe size: Child's 5½-7½	ORT29200XXS	ea.
7¼" (18.4 cm) L, 3⅝" (9.2 cm) W, X-Small, Shoe size: Women's 3-4½, Child's 10½-12½	ORT29200XS	ea.
8¾" (22.2 cm) L, 4¾" (11.1 cm) W, Small, Shoe size: Women's 5-7	ORT29200S	ea.
10½" (26.7 cm) L, 4¾" (12.1 cm) W, Medium, Shoe size: Men's 9-10, Women's 10-12	ORT29200M	ea.
12" (30.5 cm) L, 5¼" (13.3 cm) W, Large, Shoe size: Men's 12-15½	ORT29200L	ea.
13½" (34.3 cm) L, 5½" (14 cm) W, X-Large, Shoe size: Men's 16+	ORT29200XL	ea.



* Measurements based upon product size dimensions and relative shoe sizing



D Canvas Molded Rocker Cast Sandal offers a secure fit without the bulk

- Rocker sole promotes a more normal gait with heel and toe pitch
- Two canvas straps secure with hook-and-loop over the top of the foot
- Hook-and-loop heel strap for a secure fit



E Rainbow Cast Sandal is perfect for warm climates

- Multicolored elastic with hook-and-loop closures secure around the cast providing a high level of ventilation
- Thick rocker sole is made of a lightweight EVA
- Rocker sole promotes a more normal gait by incorporating both heel and toe pitch

Description	Item Number	Qty
7¼" (18.4 cm) L, 3" (7.62 cm) W, X-Small Shoe size: Child's 10½-12½	ORT290010XS	ea.
9½" (24.1 cm) L, 4¼" (10.8 cm) W, Small Shoe size: Men's 5-7½, Women's 6-8½	ORT290010S	ea.
11¼" (28.6 cm) L, 5" (12.7 cm) W, Medium Shoe size: Men's 10½-12½, Women's 11½-13	ORT290010M	ea.
12¾" (32.4 cm) L, 5½" (14 cm) W, Large Shoe size: Men's 15-16	ORT290010L	ea.
13½" (34.3 cm) L, 5½" (14 cm) W, X-Large Shoe size: Men's 16+	ORT290010XL	ea.

Description	Item Number	Qty
7¼" (18.4 cm) L, 3½" (8.9 cm) W, X-Small Shoe size: Women's 3-4½, Child's 10½-12½	ORT290040XS	ea.
9½" (24.1 cm) L, 4½" (10.5 cm) W, Small Shoe size: Men's 5-7½, Women's 5-8½	ORT290040S	ea.
10½" (26.7 cm) L, 4¾" (11.8 cm) W, Medium Shoe size: Men's 8-10, Women's 9-11	ORT290040M	ea.
11" (27.9 cm) L, 5" (12.7 cm) W, Large Shoe size: Men's 10½-12½, Women's 11½-13½	ORT290040L	ea.
12½" (31.8 cm) L, 5¾" (13.7 cm) W, X-Large Shoe size: Men's 13-16	ORT290040XL	ea.

* Measurements based upon product size dimensions and relative shoe sizing

Post-Op Shoes



- A** Vinyl Hook-and-Loop Post-Op Shoe is easy to apply
- Rigid sole prevents excessive foot motion after surgery
 - Heel and toe pitch encourages a natural gait
 - Vinyl upper resists water and wipes clean
 - Skid resistant bottom provides safe traction
 - Two hook-and-loop straps make it quick and easy to apply

Description	Item Number	Qty
Men's Shoe Size		
Small: 10-11 (Women's 11-12)	ORT30200MS	ea.
Medium: 11½-12½ (Women's 12½+)	ORT30200MM	ea.
Large: 13-14	ORT30200ML	ea.
X-Large: 14½-15½	ORT30200MXL	ea.
Women's Shoe Size		
X-Small: 3-5	ORT30200WXS	ea.
Small: Women's 6-7½	ORT30200WS	ea.
Medium: 8-9 (Men's 7-8)	ORT30200WM	ea.
Large: 9½-10½ (Men's 8½-9½)	ORT30200WL	ea.



- C** Semi-Rigid Post-Op Shoe is our most popular style
- Cushioned rocker sole with heel and toe pitch for natural gait
 - Durable nylon upper with a soft inner lining
 - Molded heel collar secures foot placement and reduces heel slipping
 - Loop-lock closure for an easy, secure fit
 - Removable tongue allows for optimum placement and comfort

Description	Item Number	Qty
Men's Shoe Size		
Small: 10-11 (Women's 11-12)	ORT30300MS	ea.
Medium: 11½-12½ (Women's 12½+)	ORT30300MM	ea.
Large: 13-14	ORT30300ML	ea.
X-Large: 14½-15½	ORT30300MXL	ea.
Women's Shoe Size		
X-Small: 3-5	ORT30300WXS	ea.
Small: 6-7½	ORT30300WS	ea.
Medium: 8-9 (Men's 7-8)	ORT30300WM	ea.
Large: 9½-10½ (Men's 8½-9½)	ORT30300WL	ea.



- B** Vinyl Lace-Up Post-Op Shoe offers protection and security
- Lace-up design provides security and a customizable size
 - Rigid sole prevents excessive foot motion after surgery
 - Heel and toe pitch encourages a natural gait
 - Vinyl upper resists water and wipes clean
 - Skid-resistant bottom provides safe traction

Description	Item Number	Qty
Men's Shoe Size		
Small: 10-11 (Women's 11-12)	ORT30100MS	ea.
Medium: 11½-12½ (Women's 12½+)	ORT30100MM	ea.
Large: 13-14	ORT30100ML	ea.
X-Large: 14½-15½	ORT30100MXL	ea.
Women's Shoe Size		
X-Small: 3-5	ORT30100WXS	ea.
Small: 6-7½	ORT30100WS	ea.
Medium: 8-9 (Men's 7-8)	ORT30100WM	ea.
Large: 9½-10½ (Men's 8½-9½)	ORT30100WL	ea.



- D** Medline's Rigid Post-Op Shoe offers protection and support during recovery
- Rigid sole reduces unwanted flexion during recovery from broken bones
 - Nylon upper with soft inner lining for patient comfort
 - Loop-lock closure for easy, secure fit

Description	Item Number	Qty
Men's Shoe Size		
Small: 10-11 (Women's 11-12)	ORT30400MS	ea.
Medium: 11½-12½ (Women's 12½+)	ORT30400MM	ea.
Large: 13-14	ORT30400ML	ea.
X-Large: 14½-15½	ORT30400MXL	ea.
Women's Shoe Size		
X-Small: 3-5	ORT30400WXS	ea.
Small: 6-7½	ORT30400WS	ea.
Medium: 8-9 (Men's 7-9)	ORT30400WM	ea.
Large: 9½-10½ (Men's 8½-9½)	ORT30400WL	ea.

* Measurements based upon relative shoe sizing