

# PROHIBIT<sup>®</sup> Series Face Masks



**An Industry First!**  
Our Exclusive Scented Face Mask  
has a Fresh Spearmint Scent.

See page 4 for details



# 5 Reasons that Make Medline's Prohibit Series® Face Masks Better than the Competition.

## 1. Incomparable Comfort and Superior Breathability

Masks that restrict air flow aren't comfortable — they're hot. Medline's PROHIBIT Series masks allow air to circulate in both directions while at the same time filtering harmful bacteria, so you stay cool, yet protected.



## 2. Bacterial Filtration Efficiency (BFE)

All PROHIBIT Series masks have bacterial filtration efficiencies exceeding 95%.\* Every mask undergoes testing by Nelson Laboratories, Inc. to ensure that you get the very best protection available.

\*Test results available upon request.

## 3. Particle Filtration Efficiency (PFE)

Surgical and Procedure masks exceed 98% PFE at 0.1  $\mu\text{m}$ .  
N95 Masks are greater than or equal to 99% PFE at 0.1  $\mu\text{m}$ .

## 4. Fluid Protection for Greater Safety

Protection from bloodborne pathogens is a serious matter. For the highest level of fluid protection, our Max X series masks pass the ASTM F1862 synthetic blood penetration at 160mm hg. These masks incorporate Drilex media, a revolutionary material prevents fluid impact penetration in either direction, yet remains very breathable.



## 5. Evaluation Program

Experience what makes Prohibit Series Face Mask better than the rest. Depending on how your account qualifies, you may be eligible to receive FREE product for evaluation. In addition, we will provide the tools and clinical support to work with staff during the evaluation. Conversion incentives may apply. Contact your Medline Sales Representative for details.





We guarantee everything we sell. Should you have any problem with a Medline product, please call us or speak to your Medline Sales Representative. We will take care of the problem right away, either by cheerfully and promptly exchanging the product or by crediting your account.

## Color Coding Index

|                              |  |
|------------------------------|--|
| Scented Masks                |  |
| Fluid Resistant Masks        |  |
| Anti-Fog Masks               |  |
| Surgical Masks               |  |
| Specialty Surgical Masks     |  |
| Isolation & Procedural Masks |  |
| N95 Particulate Respirators  |  |

**latexfree** All products in this catalog are latex-free.

## Table of Contents

|                             |       |
|-----------------------------|-------|
| Scented Masks               | 4     |
| Fluid Resistant Masks       | 4-5   |
| Anti-Fog Masks              | 6     |
| Surgical Masks              | 7     |
| Specialty Surgical Masks    | 8     |
| Isolation Masks             | 9     |
| Procedural Masks            | 9     |
| N95 Particulate Respirators | 9     |
| Cross Reference             | 10-11 |
| Mask Accessories            | 12    |



## An Industry First!

Spearmint scented mask protects against unpleasant odors in the O.R.

### Scented Procedure Mask

- Unique scent of spearmint is constant breath of fresh air.
- Fluid-resistant, outer facing provides protection.
- Soft cellulose inner facing increases inner facing comfort.
- Convenient earloops are sonically sealed to the outside of the mask to virtually eliminate possible facial irritation and ensure a secure closure.
- Full-width nosepiece guarantees proper fit.

| Item #   | Description                                 | Pkg           |
|----------|---|---------------|
| NON27605 | Scented Procedure Mask with Earloops, Green | 50/bx, 300/cs |

BFE > 97%  
PFE > 98%



### Scented Surgical Mask

- Unique scent of spearmint is constant breath of fresh air.
- Fluid-resistant, outer facing provides protection.
- Soft cellulose inner facing increases inner facing comfort.
- Full-width nosepiece guarantees proper fit.

| Item #   | Description                            | Pkg           |
|----------|--|---------------|
| NON27610 | Scented Surgical Mask with Ties, Green | 50/bx, 300/cs |

BFE > 97%  
PFE > 98%



## Max X Fluid Protection

Our Highest Level of Fluid Protection Ever!

### Max X Fluid Protection Masks w/Eyeshield

- Non-glare, fog-free, wraparound eyeshield assures maximum eye and facial protection.
- Unique Drilex filter provides uncompromising fluid resistance.
- Cellulose inner and outer facing offers protection against fluid and inner facing comfort.
- Anti-fog strip of comfortable, soft foam prevents eyeglasses from fogging.

| Item #     | Description                             | Pkg           |
|------------|---|---------------|
| NON27710   | Max X w/Eyeshield and Ties, Pattern     | 25/bx, 100/cs |
| NON27710EL | Max X w/Eyeshield and Earloops, Pattern | 25/bx, 100/cs |

BFE > 99%  
PFE > 99%

ASTM F1862 160mm hg  
Fluid Resistance



### Max X Fluid Protection Masks

- Cellulose inner and outer facing offers protection against fluid and inner facing comfort.
- Unique Drilex filter provides uncompromising fluid resistance.
- Anti-fog strip of comfortable, soft foam prevents eyeglasses from fogging, thus ensuring optical clarity.
- Full-width nosepiece guarantees proper fit.

| Item #     | Description                  | Pkg           |
|------------|------------------------------|---------------|
| NON27712   | Max X with Ties, Pattern     | 50/bx, 300/cs |
| NON27712EL | Max X with Earloops, Pattern | 50/bx, 300/cs |

BFE > 99%  
PFE > 99%

ASTM F1862 160mm hg  
Fluid Resistance



## Fluid Resistant Masks

Great fluid protection while filtering out bacteria

BFE > 99%  
PFE > 99%

ASTM F1862 120mm hg  
Fluid Resistance

### Fluid-Resistant Masks with Eyeshield

- Non-glare, fog-free, wraparound eyeshield assures maximum eye and facial protection.
- Polypropylene outer facing for added protection against fluids, plus soft inner cellulose facing for comfort.
- Anti-fog strip of comfortable, soft foam prevents eyeglasses from fogging.

| Item #     | Description   | Pkg           |
|------------|---|---------------|
| NON27410   | Fluid Resistant with Eyeshield and Ties, Purple     | 25/bx, 100/cs |
| NON27410EL | Fluid Resistant with Eyeshield and Earloops, Purple | 25/bx, 100/cs |



BFE > 99%  
PFE > 99%

ASTM F1862 120mm hg  
Fluid Resistance

### Fluid-Resistant Masks

- Polypropylene outer facing for added protection against fluids, plus soft inner cellulose facing for comfort.
- Anti-fog strip of comfortable, soft foam prevents eyeglasses from fogging, thus ensuring optical clarity.
- Full-width nosepiece guarantees proper fit.

| Item #     | Description                           | Pkg           |
|------------|---------------------------------------|---------------|
| NON27412   | Fluid Resistant with Ties, Purple     | 50/bx, 300/cs |
| NON27412EL | Fluid Resistant with Earloops, Purple | 50/bx, 300/cs |



## Prohibit X-Tra™

Offers More Fluid Protection than Traditional Surgical Masks

BFE > 98%  
PFE > 98%

ASTM F1862 80mm hg  
Fluid Resistance

### Prohibit X-Tra™ Fluid Protection Masks w/Eyeshield

- Non-glare, fog-free, wraparound eyeshield assures maximum eye and facial protection.
- Cellulose inner and outer facing offers protection against fluid and inner facing comfort.
- Anti-fog strip of comfortable, soft foam prevents eyeglasses from fogging.

| Item #   | Description                                    | Pkg           |
|----------|--|---------------|
| NON27405 | X-Tra Protection with Eyeshield and Ties, Blue | 25/bx, 100/cs |



BFE > 98%  
PFE > 98%

ASTM F1862 80mm hg  
Fluid Resistance

### Prohibit X-Tra™ Fluid-Protection Masks

- Cellulose inner and outer facing offers protection against fluid and inner facing comfort.
- Anti-fog strip of comfortable, soft foam prevents eyeglasses from fogging, thus ensuring optical clarity.
- Full-width nosepiece guarantees proper fit.

| Item #     | Description                          | Pkg           |
|------------|--------------------------------------|---------------|
| NON27408   | X-Tra Protection with Ties, Blue     | 50/bx, 300/cs |
| NON27408EL | X-Tra Protection with Earloops, Blue | 50/bx, 300/cs |



## Anti-Fog Mask Styles

With four different ways to protect you from "fogging up" Medline gives you more choices.



### Anti-Fog Masks with Adhesive Tape

- Two styles to choose from: standard adhesive or aggressive adhesive.
- Non-irritating, adhesive strip forms a tight seal to prevent eyeglass fogging and maintain optical clarity.
- Fluid-resistant facing provides excellent protection.
- Full-width nosepiece guarantees proper fit.

BFE > 97%  
PFE > 98%

| Item #    | Description                          | Pkg           |
|-----------|--------------------------------------|---------------|
| NON27371A | Adhesive with Ties, Green            | 50/bx, 300/cs |
| NON27361A | Aggressive Adhesive with Ties, Green | 50/bx, 300/cs |

### Anti-Fog Mask with Film Strip

- A film strip is designed to conform to the face to help seal in air, reducing eyeglass and microscope fogging.
- Cellulose facing for cool comfort.
- Full-width nosepiece guarantees proper fit.

BFE > 97%  
PFE > 98%

| Item #    | Description                                   | Pkg           |
|-----------|---|---------------|
| NON27379A | Anti-Fog Mask with Film Strip and Ties, Green | 50/bx, 300/cs |

### Anti-Fog Mask with Foam Strip

- Thick, soft foam strip keeps condensation from escaping out the top of the mask thus preventing eyeglass fogging and ensuring optical clarity.
- Cellulose facing for cool comfort.
- Full-width nosepiece guarantees proper fit.

BFE > 97%  
PFE > 98%

| Item #    | Description                                  | Pkg           |
|-----------|--|---------------|
| NON27378A | Anti-Fog Mask with Foam Strip and Ties, Blue | 50/bx, 300/cs |



### Anti-Fog Chamber-Style Mask with Foam Strip

- Large breathing chamber keeps mask away from face for greater comfort and easier breathing; also allows unmuffled speech.
- Thick, foam strip keeps condensation from escaping out the top of the mask, thus preventing eyeglass fogging and ensuring optical clarity.
- Polypropylene outer facing for extra fluid protection along with soft cellulose inner facing for cool comfort.

BFE > 97%  
PFE > 98%

| Item #    | Description                                     | Pkg           |
|-----------|---|---------------|
| NON27373A | Chamber Style Mask w/Foam Strip and Ties, Green | 50/bx, 300/cs |



## Surgical Masks

Our Surgical Masks Give you Peace of Mind in the O.R.

### Standard Surgical Mask for your every day needs.

- NON27376 and NON27600 features a cellulose inner and outer facing keeps the wearer cool and comfortable.
- NON27377 and NON27420 feature a soft, fluid resistant polypropylene inner and outer facing assures comfort on the inside and protection on the outside.
- Full-width nosepiece guarantees proper fit.

| Item #   | Description                              | Pkg           |
|----------|--|---------------|
| NON27376 | Standard Surgical Mask and Ties, Blue    | 50/bx, 300/cs |
| NON27600 | Horizontal Tie Surgical Mask, Blue       | 50/bx, 300/cs |
| NON27377 | Surgical Mask and Ties, Blue             | 50/bx, 300/cs |
| NON27420 | Surgical Mask with Shield and Ties, Blue | 25/bx, 200/cs |

BFE > 97%  
PFE > 98%



### Medsoft™ Surgical Mask is our Softest.

- Our softest mask available.
- Fluid-resistant, outer facing provides protection while the soft inner facing gives the wearer the most comfort against the face.
- Full-width nosepiece guarantees proper fit.

| Item #   | Description                        | Pkg           |
|----------|------------------------------------|---------------|
| NON27400 | Soft Surgical Mask and Ties, White | 50/bx, 300/cs |

BFE > 95%  
PFE > 98%



### MEDLITE™ Surgical Mask feels so light you'll forget your wearing it.

- Our most comfortable and cool mask, ideal for long procedures.
- Fluid-resistant, outer facing provides protection.
- Soft cellulose inner facing keeps the wearer cool and comfortable.
- Full-width nosepiece guarantees proper fit.

| Item #   | Description                          | Pkg           |
|----------|--------------------------------------|---------------|
| NON27402 | MEDLITE Surgical Mask and Ties, Blue | 50/bx, 300/cs |

BFE > 97%  
PFE > 98%



### Cone-Style Mask Enables Easy Donning

- Fluid-resistant, cup-shaped mask is held away from the nose and mouth for comfort and protection.
- Cone-style enables easy donning.
- Malleable nosepiece allows a custom fit.
- Designed with comfortable, rounded edges.

| Item #   | Description                         | Pkg           |
|----------|-------------------------------------|---------------|
| NON27381 | Cone Style Mask with Headband, Blue | 50/bx, 300/cs |

BFE > 99%  
PFE > N/A



### Chamber-Style Mask Features a Large Breathing Chamber for Greater Comfort

BFE > 97%  
PFE > 98%

- Large breathing chamber keeps mask away from face for greater comfort and easier breathing; also allows unmuffled speech.
- An ideal mask for beginning health care professionals as well as veterans.
- Cellulose facing for cool comfort.
- Full-width nosepiece guarantees proper fit.

| Item #   | Description                        | Pkg           |
|----------|------------------------------------|---------------|
| NON27382 | Chamber Style Mask and Ties, Green | 50/bx, 300/cs |



### Chamber-Style Mask with Eyeshield Assures Maximum Eye and Facial Protection

- Large breathing chamber keeps mask away from face for greater comfort
- Non-glare, fog-free, wraparound eyeshield assures maximum eye and facial protection.
- Polypropylene outer facing for added protection against fluids, plus soft inner cellulose facing for comfort.
- Anti-fog strip of comfortable foam prevents eyeglasses from fogging.

| Item #   | Description                                    | Pkg           |
|----------|--|---------------|
| NON27411 | Chamber Style Mask w/Eyeshield and Ties, Green | 25/bx, 100/cs |

BFE > 97%  
PFE > 98%

### Hypoallergenic Mask is Ideal for Sensitive Skin

- Contains no dyes or colorants that could irritate sensitive skin.
- Comfortable cellulose outer and inner facing
- Full-width nosepiece guarantees proper fit.

| Item #   | Description                          | Pkg           |
|----------|--------------------------------------|---------------|
| NON27385 | Hypoallergenic Mask with Ties, White | 50/bx, 300/cs |

BFE > 97%  
PFE > 98%



### Laser Mask Protects Against Laser Plume Particles

- Special filter blocks laser plume particles as small as 0.1 micron.
- Fluid-resistant outer facing combined with comfortable, soft cellulose inner facing.
- Full-width nosepiece guarantees proper fit.

| Item #   | Description                | Pkg           |
|----------|----------------------------|---------------|
| NON27374 | Laser Mask with Ties, Gray | 50/bx, 300/cs |

BFE > 99%  
PFE > 99%



### Did you know?

Medline also offers a full line of disposable protective apparel that complements our fine face masks. Ask your Medline sales representative for details.



### Isolation Mask

- Fluid-resistant outer facing provides protection.
- Soft cellulose inner facing increases comfort on the inside.
- Full-width nose-piece guarantees proper fit.

| Item #   | Description                          | Pkg           |
|----------|--------------------------------------|---------------|
| NON27120 | Isolation Mask with Earloops, Yellow | 50/bx, 300/cs |

BFE > 97%  
PFE > 98%



### Basic Isolation Mask

- NON27110 has a tissue inner and outer facing making mask light and comfortable.
- NON27122 features a soft, fluid resistant polypropylene inner and outer facing assures comfort on the inside and protection on the outside.
- Full-width nose-piece guarantees proper fit.

| Item #   | Description                                | Pkg           |
|----------|--|---------------|
| NON27110 | Basic Isolation Mask with Earloops, Yellow | 50/bx, 300/cs |
| NON27122 | Isolation Mask with Earloops, Yellow       | 50/bx, 300/cs |

BFE > 99%  
PFE > N/A

### Standard Procedure Mask

- NON27375 features a Soft cellulose inner facing increasing comfort on the inside.
- NON27420EL features a soft, fluid resistant polypropylene inner and outer facing assures comfort on the inside and protection on the outside.
- Full-width nose-piece guarantees proper fit.

| Item #     | Description                                   | Pkg           |
|------------|---|---------------|
| NON27375   | Standard Procedure Mask with Earloops, Blue   | 50/bx, 300/cs |
| NON27420EL | Procedure Mask with Shield and Earloops, Blue | 25/bx, 200/cs |

BFE > 95%  
PFE > 98%



### Basic Procedure Mask

- NON27365 has a tissue inner and outer facing making mask light and comfortable.
- NON27378 features a Soft, fluid resistant polypropylene inner and outer facing assures comfort on the inside and protection on the outside.
- Full-width nose-piece guarantees proper fit.

| Item #   | Description                              | Pkg           |
|----------|--|---------------|
| NON27365 | Basic Procedure Mask with Earloops, Blue | 50/bx, 300/cs |
| NON27378 | Basic Procedure Mask with Earloops, Blue | 50/bx, 300/cs |

BFE > 99%  
PFE > N/A

## N95 Approved Particulate Respirators Meet NIOSH and CDC Guidelines for TB Exposure

### Adjustable, Flat Fold and Cone Styles

- N95 approved! Meets NIOSH and CDC guidelines for TB exposure.
- Positive Facial Lock® and malleable nose-piece allow a custom fit for any face, so one size mask fits all. (NON27501 only)
- Micro filter medium protects against micron-size particles.
- Double band holds mask firmly in place for maximum protection.
- All respirators must be selected, fitted and maintained in accordance with OSHA and other applicable regulations.

| Item #   | Description  | Pkg           |
|----------|--|---------------|
| NON27501 | Adjustable Particulate Respirator, Green/White Pattern | 35/bx, 210/cs |
| NON24506 | Cone Style Particulate Respirator, Large/Medium, White | 20/bx, 240/cs |
| NON24507 | Cone Style Particulate Respirator, Small, White        | 20/bx, 240/cs |
| NON24505 | Flat Fold Respirator, Dual Headband, White             | 20/bx, 200/cs |



## Prohibit Series Competitive Cross Reference

| Description                              | Medline Item # | Medline Packaging | Kimberly-Clark  |
|--|----------------|-------------------|---|
| Isolation mask                           | NON27110       | 50/bx, 300/cs     | 47117   |
|  | NON27122       | 50/bx, 300/cs     |   |
|  | NON27120       | 50/bx, 300/cs     |   |
| Procedure mask                           | NON27365       | 50/bx, 300/cs     | 47080<br>62356, 62357   |
|  | NON27375       | 50/bx, 300/cs     |   |
|  | NON27378       | 50/bx, 300/cs     |   |
|  | NON27420EL     | 25/bx, 100/cs     |   |
|  | NON27605       | 50/bx, 300/cs     |   |
| Standard surgical                        | NON27376       | 50/bx, 300/cs     | 48201, 48210<br>47500<br>48100, 48105<br>48390  |
|  | NON27600       | 50/bx, 300/cs     |   |
|  | NON27402       | 50/bx, 300/cs     |   |
|  | NON27400       | 50/bx, 300/cs     |   |
|  | NON27377       | 50/bx, 300/cs     |   |
|  | NON27420       | 25/bx, 100/cs     |   |
|  | NON27610       | 50/bx, 300/cs     |   |
| Anti-Fog                                 | NON27361A      | 50/bx, 300/cs     | 49215<br>49235<br>49310, 49214  |
|  | NON27371A      | 50/bx, 300/cs     |   |
|  | NON27379A      | 50/bx, 300/cs     |   |
|  | NON27378A      | 50/bx, 300/cs     |   |
| Chamber style                            | NON27382       | 50/bx, 300/cs     | 48220, 37525, 48227<br>49216<br>48228   |
|  | NON27373A      | 50/bx, 300/cs     |   |
|  | NON27411       | 25/bx, 100/cs     |   |
| Hypoallergenic                           | NON27385       | 50/bx, 300/cs     | 59928   |
| Cone                                     | NON27381       | 50/bx, 300/cs     | 152   |
| Laser                                    | NON27374       | 50/bx, 300/cs     | 47625, 47650  |
| Fluid protection                         | NON27710       | 25/bx, 100/cs     | 48247, 48237, 48248<br>47137, 47147<br>48207, 48208, 48297<br>47107, 47297<br>62114<br>62116, 47298<br>62113<br>62115 |
|  | NON27710EL     | 25/bx, 100/cs     |   |
|  | NON27712       | 50/bx, 300/cs     |   |
|  | NON27712EL     | 50/bx, 300/cs     |   |
|  | NON27410       | 25/bx, 100/cs     |   |
|  | NON27410EL     | 25/bx, 100/cs     |   |
|  | NON27412       | 50/bx, 300/cs     |   |
|  | NON27412EL     | 50/bx, 300/cs     |   |
|  | NON27405       | 25/bx, 100/cs     |   |
|  | NON27408       | 50/bx, 300/cs     |   |
|  | NON27408EL     | 50/bx, 300/cs     |   |
| N95 Particulate respirator & accessories | NON24506       | 20/bx, 240/cs     | 46727, 46767, 62126<br>46827, 46867, 62355<br>47950   |
|  | NON24507       | 20/bx, 240/cs     |   |
|  | NON24505       | 50/bx, 300/cs     |   |
|  | NON27501       | 35/bx, 210/cs     |   |
|  | NON24510       | each              |   |
| Face shields and protective eyewear      | NONFS300       | 24/bx, 96/cs      | 41204<br>41205<br>SV250L<br>SV100F  |
|  | NONFS400       | 24/bx, 96/cs      |   |
|  | NONLENS        | 250/bx            |   |
|  | NONFRAME       | 100/bx            |   |
|  | NON24774       | 12/bx             |   |
|  | NON24777A      | 20/bx             |   |
|  | NON24777S      | 20/bx             |   |
|  | NON24775       | 36/cs             |   |
| NON24776                                 | 36/cs          |                   |   |

## Prohibit Series Competitive Cross Reference

| Medline Item #   | 3M  | Cardinal  | Precept   |
|--|---|---|---|
| NON27110<br>NON27122<br>NON27120   |   | AT70061,AT70021<br>AT7004                                   | 15100, 15120<br>15110   |
| NON27365<br>NON27375<br>NON27378<br>NON27420EL<br>NON27605   | 1826<br>1820  | AT71061,AT71021<br><br>AT7511<br>AT7511-WE                  | 15101<br><br>15111  |
| NON27376<br>NON27600<br>NON27402<br>NON27400<br>NON27377<br>NON27420<br>NON27610   | 1818<br>1816  | AT7505<br>AT71039<br>AT71035<br><br><br>AT71135,AT71335     | 15201, 65-3110<br>15250<br>15200, 65-3010<br><br>15600, 65-3100 |
| NON27361A<br>NON27371A<br>NON27379A<br>NON27378A   | 1831<br>1834<br>1832                                    | AT72835<br>AT752005<br>AT71235,AT72235                      | 15212, 65-3321, 65-3322<br>15211<br>15210, 65-3320              |
| NON27382<br>NON27373A<br>NON27411  | 1836, 1838<br>1838FSG                                   | AT51035<br>AT54535<br>AT51135,AT54635                       | 15220, 65-3124<br><br>65-3346                                   |
| NON27385   | 1833  | AT73035   | 15215   |
| NON27381   | 1800ml  | AT7509  | 65-3333   |
| NON27374   |   | AT75435   | 65-3310, 65-3311  |
| NON27710<br>NON27710EL<br>NON27712<br>NON27712EL<br>NON27410<br>NON27410EL<br>NON27412<br>NON27412EL<br>NON27405<br>NON27408<br>NON27408EL | 1835FS<br>1840FS<br>1835<br>1840<br><br><br><br>1818FSG | AT74635<br>AT74631<br>AT74535<br>AT74531                    | 15330<br>15310, 15510<br>15320<br>15300, 15301, 15500           |
| NON24506<br>NON24507<br>NON24505<br>NON27501<br>NON24510   | 1860<br>1860S<br>1870<br><br>FT10                       | N95-ML<br>N95-S   | 65-3395<br>65-3395s<br><br>65-65000                             |
| NONFS300<br>NONFS400<br>NONLENSE<br>NONFRAME<br>NON24774<br>NON24777A<br>NON24777S<br>NON24775<br>NON24776                                 | <br><br><br><br>1710<br>1611<br><br>1621                | FISHIELD50<br>HISHIELD50<br>ESHIELD50<br>BFRAME30, WFRAME30 | 15995<br>15990  |

## Face Mask Accessories To Complement Our Face Masks

### Qualitative Fit Test Kit

- To be used in conjunction with our N95 particulate respirator masks (NON27501, NON24505, NON24506 and NON24507).
- Fit testing is required prior to the initial use of a respirator.
- Saccharin based fit test kit.
- Test is easy to administer and instructions are found within the kit.
- One kit tests approximately 100 staff members.

| Item #    | Description                   | Pkg      |
|-----------|-------------------------------|----------|
| NON24510  | Qualitative Fit Test Kit      | each     |
| NON24510A | Fit Test Sensitivity Solution | 6 btl/cs |
| NON24510B | Fit Test Actual Solution      | 6 btl/cs |
| NON24510C | Fit Test Bitter Solution      | 6 btl/cs |
| NON24510D | Fit Test Hood,                | 10/cs    |
| NON24510E | Replacement Nebulizers        | 3/cs     |



### Disposable Face Shields

- Provides the user with maximum eye and facial protection.
- Can be used in conjunction with our Prohibit Face Masks.
- Shields are made of durable 6 ml thick material and are anti-fog to ensure optical clarity.
- Both shields have a latex-free band and comfortable, contoured foam material that rests against the forehead.

| Item #   | Description                       | Pkg          |
|----------|-----------------------------------|--------------|
| NONFS300 | Disposable Full Shield, 7.5" long | 24/bx, 96/cs |
| NONFS400 | Disposable 3/4 Shield, 6.5" long  | 24/bx, 96/cs |



### Fluid Protective Goggles

- Soft, vinyl-frame goggles fit easily over prescription eyewear.
- Four air vents provide additional comfort.
- Available with specially-treated, anti-fog lenses for clear, unobstructed viewing.
- Goggles meet ANSI standards and help your facility meet OSHA requirements for hazardous liquids.

| Item #   | Description       | Pkg   |
|----------|-------------------|-------|
| NON24775 | Goggles, fog free | 36/cs |
| NON24776 | Goggles, standard | 36/cs |



### Safety Glasses

- Safety glasses provide economical and practical protection.
- Built-in brow bar and side shields, plus wide view, optical lenses.
- Comfortable fit.
- Large size glasses can be worn over prescription eyewear.

| Item #    | Description                                     | Pkg    |
|-----------|---|--------|
| NON24777A | Safety glasses, visitor, clear, large           | 20/bx  |
| NON24777S | Safety glasses, visitor, clear, small           | 20/bx  |
| NON24774  | Safety glasses with black strip                 | 12/bx  |
| NONLENS   | Disposable lens for use w/NONFRAME, clear       | 250/bx |
| NONFRAME  | Reusable frame for use NONLENS, assorted colors | 100/bx |



### Face Mask Dispensers

| Item #   | Description                        | Pkg  |
|----------|------------------------------------|------|
| NONFP022 | Dispenser for surgical face masks  | 6/cs |
| NONFP021 | Dispenser for procedure face masks | 6/cs |



1-800-MEDLINE (633-5463) [www.medline.com](http://www.medline.com)

©2008 Medline Industries, Inc. Prohibit and Medline are registered trademarks of Medline Industries, Inc. One Medline Place, Mundelein, IL 60060. Positive Facial Lock is a registered trademark of Alpha ProTech, Inc. We reserve the right to correct any errors that may occur within this brochure.

Medline Canada 1-800-396-6996

[www.medline.ca](http://www.medline.ca) [canada@medline.com](mailto:canada@medline.com)

Not all products are currently available for sale in Canada. Please contact your Account Representative for additional information. MKT208031/LIT168/15M/ZPN6