

MILTEX HI-LIGHT

TIP-IT® INSTRUMENT GUARDS



Miltex® Tip-It® Instrument Guards protect instrument tips and edges from damage during handling and sterilization. Miltex now offers a newly expanded line of Tip-It® Instrument Guards with three additional styles: Clear Vented, Tinted Vented, and Solid Vented in addition to the Solid Non-Vented style. Miltex® Tip-It® Instrument Guards feature:

FEATURES

- PVC medical grade plastic for use in medical applications
- Latex-free material that is non-allergenic
- Large holes (vented styles) to reduce vacuum and moisture build up for easier removal from instrument
- Steam and gas permeable to allow for effective sterilization of instruments
- Four styles in nine colors and sizes for a variety of choices
- Non-sterile disposable, 50/pkg, including assortment pack

Part numbers and descriptions on reverse.

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.

Miltex
Redefining Excellence

589 Davies Drive
York, PA 17402
phone 717 840-9335
toll free 866 854-8300
fax 866 854-8400
www.miltex.com

0107/5/C/3

SURG-TIPIT-HI





REF Clear Vented	REF Tinted Vented	REF Solid Vented	REF Solid Non-Vented	Dimension	Size Code	Color
3-2501C	3-2501T	3-2501V	3-2501	1/16" X 3/4" (1.6mm X 19mm)	1	White*
3-2502C	3-2502T	3-2502V	3-2502	5/64" X 3/4" (2mm X 19mm)	2	Blue*
3-2503C	3-2503T	3-2503V	3-2503	7/64" X 3/4" (2.8mm X 19mm)	3	Green*
3-2504C	3-2504T	3-2504V	3-2504	1/8" X 1" (3.2mm X 25.4mm)	4	Red*
3-2505C	3-2505T	3-2505V	3-2505	3/16" X 1" (4.8mm X 25.4mm)	5	Yellow*
3-2506C	3-2506T	3-2506V	3-2506	13/32" X 3/4" (10.3mm X 19mm)	6	Maroon*
3-2507C	3-2507T	3-2507V	3-2507	1/16" X 3/8" X 1" (1.6mm X 9.5mm X 25.4mm)	7	Orange*
3-2508C	3-2508T	3-2508V	3-2508	5/64" X 5/8" X 1" (2mm X 16mm X 25.4mm)	8	Brown*
3-2509C	3-2509T	3-2509V	3-2509	1/8" X 1" X 1" (3.2mm X 25.4mm X 25.4mm)	9	Black*
3-2510C	3-2510T	3-2510V	3-2510	Assorted, equal number of sizes plus 5		

* Color does not apply to clear vented pattern