


Premium Vinyl Exam Gloves

Advanced Technology, for all-around better performance!

Three Generations:

- **First:** standard vinyl material with PVC resin and plasticizers
- **Second:** stretch vinyl material with softer PVC resins for better fit, feel and performance.
- **Third:** patented stretch vinyl material with ultra-high-molecular-weight resins and softening and strengthening agents for extraordinary tear resistance, comfort and barrier protection.



[Click to Enlarge](#)

Generation Pink Vinyl



- The best performing vinyl glove on the market
- Made with patented 3rd generation vinyl technology
- **\$1.00/case of gloves sold will be donated to the National Breast Cancer Foundation**
- 6.7 mil finger thickness; 5.5 mil palm thickness

Aloetouch® Ultra IC Stretch Vinyl



- Inside coated with pure organic aloe to moisturize skin
- Made with 2nd-generation vinyl technology
- 6.7mil finger thickness; 5.5mil palm thickness

Aloetouch® 3G Stretch Vinyl



- Inside coated with pure organic aloe to moisturize skin
- The best performing vinyl glove on the market
- Made with 3rd-generation vinyl technology
- 6.7 mil finger thickness; 5.5 mil palm thickness

CURAD® 3G Stretch Vinyl



- The brand you know and trust
- Improved material offers remarkable softness, strength and elasticity
- Made with 3rd-generation vinyl technology
- 6.7 mil finger thickness; 5.5 mil palm thickness

CURAD® Stretch Vinyl



- The brand you know and trust
- Improved material offers a good fit and tactile sensitivity
- Made with 2nd-generation vinyl technology
- 5.9 mil finger thickness; 5.1 mil palm thickness

Ultra Stretch Vinyl



- Advanced formulation synthetic vinyl stretches to fit your hand
- Improved material offers a good fit and tactile sensitivity
- Made with 2nd-generation vinyl technology
- 6.3 mil finger thickness; 5.1 mil palm thickness

Accutouch Vinyl



- Cost-effective alternative, perfect for general patient care
- Made with 1st-generation vinyl technology
- 5.9 mil finger thickness; 5.1 mil palm thickness

MediGuard® Vinyl



- Cost-effective alternative, perfect for general patient care
- Made with 1st-generation vinyl technology
- 5.9 mil finger thickness; 4.7 mil palm thickness

Ultra 12" Stretch Vinyl



- Advanced formulation synthetic vinyl stretches to fit your hand
- Improved material offers remarkable softness, strength and elasticity
- Made with 2nd-generation vinyl technology
- 2.8 mil finger thickness; 3.1 mil palm thickness

Sterile Vinyl Exam Glove Options



- Advanced formulation for better tactile sensitivity and fit
- Plastic packaging resists moisture and minimizes fiber tears
- Small-sized packaging is easy to open and takes up less storage space
- Smooth — not textured
- 5.9 mil finger thickness; 5.1 mil palm thickness
- Available in singles or pairs