

## TRANSEPTIC® CLEANING SOLUTION

Multi-purpose spray for professional use in cleansing ultrasound transducer/probe surfaces between procedures or as a cleaning step prior to sterilization. Will not cause transducer membrane swelling.

### Features

- Convenient clear spray dispenser
- Does not contain phenol or glutaraldehyde
- Accepted for use by manufacturers of ultrasound equipment

Product #	Package Size
-----------	--------------

09-25	250 ml (8.5 fl. oz.) clear spray bottle, 12 per box (Ground shipment only)
-------	---