



**ULTRARAY**  
MEDICAL PRODUCTS INC.

**RADIATION PROTECTION & ACCESSORIES**

[HOME PAGE](#)

TEL: (905) 338-6857  
FAX: (866) 338-8990  
TOLL FREE: (877) 338-6857

 **PROUD TO BE CANADIAN MADE** 

Since its inception in 1983 UltraRay® has continued to be an industry leader in the field of radiation protection. We manufacture Canadian made products using the latest technology while maintaining the highest safety standards. Our products are manufactured using the strictest quality control procedures in accordance with Health Canada Medical Devices Regulations. We are a Class 1 Certified manufacturer and distributor operating under our Medical Establishment License 2312.

Our distributor network of 50 sales representatives across Canada provide on-site services such as samples, measurement and needs analysis consultation to ensure that clients receive the products best suited to their needs. Client input is utilized in the ongoing improvement & development of our products.

All our products carry a 2 year Warranty and 100% satisfaction guarantee.

Copyright © 2005-2010 UltraRay Medical Products Inc., All Rights Reserved

## CLEAR Pb® SHIELDING AND MOBILE BARRIERS

Mobile X-Ray Barriers shield medical personnel from secondary radiation in Radiology, Diagnostic Imaging, Nuclear Medicine, Cath Lab, O.R. and special procedure rooms.

### Features:

- Made with Clear Pb® Lead Acrylic Windows that are Shatter Resistant
- Variety of Sizes Available for Shielding 1, 2 or More People
- 8 Colours Available ( Black, Grey, Beige, White, Blue, Green, Yellow, Red )
- Heavy Duty Powder Coated Steel Construction
- Available in 0.5mm / 0.8mm / 1.0mm / 2.0mm Lead Equivalencies
- WE CUSTOM BUILD TO YOUR SPECIFICATIONS
- LIMITED LIFETIME WARRANTY



### CLEAR-Pb® - Mobile X-Ray Barrier specifications

	Regular UR 603	Jumbo UR 604	Deluxe UR 605      UR 606		Ultra UR 607      UR 608		Nuclear Medicine UR 609
<b>Overall size</b>	75 (h) x 30 in (w)	75 (h) x 30 in (w)	75 (h) x 48 in (w)	75 (h) x 48 in (w)	75 (h) x 72 in (w)	75 (h) x 72 in (w)	61 (h) x 36 in (w)
<b>Window pane lead equivalent</b>	24 (h) x 30 in (w) 0.5 mm	60 (h) x 30 in (w) 0.5 mm	36 (h) x 48 in (w) 0.5 mm	36 (h) x 48 in (w) 1.0 mm	36 (h) x 72 in (w) 0.5 mm	36 (h) x 72 in (w) 1.0 mm	24 (h) x 36 in (w) 0.8 mm
<b>Opaque panel lead equivalent</b>	48 (h) x 30 in (w) 0.06 in (1.5 mm)	12 (h) x 30 in (w) —	36 (h) x 48 in (w) 0.06 in (1.5 mm)	36 (h) x 48 in (w) 0.06 in (1.5 mm)	36 (h) x 72 in (w) 0.06 in (1.5 mm)	36 (h) x 72 in (w) 0.06 in (1.5 mm)	34 (h) x 36 in (w) 0.06 in (1.5 mm)
<b>Shipping weight</b>	110 lb	113 lb	230 lb	312 lb	334 lb	409 lb	170 lb

## CLEAR-Pb® Lead-Plastic X-Ray Shielding Window Panes

Any size window pane can be used in a thru-wall installation. We offer a complete line of Lead-Lined Steel Window Frames for use in thru-wall applications.

Window size (inches)	Lead equivalence 0.3 mm (thickness 8 mm)	Lead equivalence 0.5 mm (thickness 12 mm)	Lead equivalence 0.8 mm (thickness 18 mm)	Lead equivalence 1.0 mm (thickness 22 mm)	Lead equivalence 1.5 mm (thickness 35 mm)	Lead equivalence 2.0 mm (thickness 46 mm)
	Model	Model	Model	Model	Model	Model
<b>12 x 12</b>		<b>56-741</b>	<b>56-751</b>	<b>56-761</b>	<b>56-771</b>	<b>56-791</b>
12 x 24		56-740	56-750	56-760	56-770	56-790
<b>18 x 24</b>		<b>56-708</b>	<b>56-718</b>	<b>56-728</b>	<b>56-738</b>	
18 x 48		56-702	56-712	56-722	56-732	
<b>24 x 24</b>	<b>56-786</b>	<b>56-705</b>	<b>56-715</b>	<b>56-725</b>	<b>56-735</b>	<b>56-796</b>
24 x 30		56-742	56-752	56-762	56-772	
<b>24 x 36</b>		<b>56-704</b>	<b>56-714</b>	<b>56-724</b>	<b>56-734</b>	
24 x 48	56-787	56-747	56-757	56-767	56-777	56-797
<b>32 x 40</b>		<b>56-701</b>	<b>56-711</b>	<b>56-721</b>	<b>56-731</b>	
36 x 48	56-788	56-703	56-713	56-723	56-733	56-798
<b>36 x 60</b>		<b>56-706</b>	<b>56-716</b>	<b>56-726</b>	<b>56-736</b>	
36 x 72		56-707	56-717	56-727	56-737	
<b>36 x 84</b>	<b>56-780</b>	<b>56-709</b>	<b>56-719</b>	<b>56-729</b>	<b>56-739</b>	<b>56-799</b>
48 x 48		56-746	56-756	56-766	56-776	56-792
<b>48 x 60</b>		<b>56-745</b>	<b>56-755</b>	<b>56-765</b>	<b>56-775</b>	
48 x 72	56-789	56-700	56-710	56-720	56-730	
<b>48 x 84</b>		<b>56-744</b>	<b>56-754</b>	<b>56-764</b>	<b>56-774</b>	
48 x 96		56-743	56-753	56-763	56-773	56-793
<b>72 x 96</b>		<b>56-748</b>	<b>56-758</b>	<b>56-768</b>	<b>56-778</b>	

**NOTE: WE CAN CUSTOM BUILD MOBILE BARRIERS FOR ALMOST ALL OF THE ABOVE PANES**

## MOBILE SHIELDS



Excellent for patient protection during diagnostic procedures such as erect chest and cervical spine. Includes a durable nylon covered 0.5mm shield. Height adjustable 36" – 60". Rolls on a sturdy star-leg aluminum base for extra stability and easy movement.

MSB-2 MOBILE SOFT BARRIER

## ADJUSTABLE HEIGHT SHIELD



Designed to protect the patient against scattered radiation. Ideal for chest, skull, and spine radiography. Made of lead vinyl with .5mm lead protection and attractively covered on both sides.

26" (65cm) wide and adjustable with a minimum height of 38" and a maximum height of 58". Mounted on welded durable chrome frame and rolls easily on 4 casters.


Specifications: 58" H (open) 38" H (closed) x 26" W x 18" D (at base)

MSB-1 ADJUSTABLE HEIGHT SHIELD

## FLUOROSCOPY SCREENS / TABLE SHIELDS

UltraRay manufactures low cost replacement fluoroscopy screens and shields for most standard tables.

We can custom make all sizes and types with or without hardware at well below manufacturer prices. Send us your shields for recovery or contact us to discuss your custom needs.





## CURTAINS

Curtains are ideal for areas where secondary or low level radiation is present. Optional hardware or grommets can be added.

- Doorway Curtains
- Room Divider Curtains
- 0.5mm or 1.0mm Lead Equivalency Available
- Custom Made to Your Needs



## BLANKETS

Our lead blankets are made of lead wool or lead vinyl depending on the application. All blankets are custom made to your specifications.

- Box Blankets
- Equipment Blankets
- 0.5mm, 1.0mm, 1.5mm and 2.0mm LE
- Available in 30 Colours